Le 26 juin 2024

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
Référé par:

Non validé

# Fluke Biomedical Ansur Test and Inspection Procedure

Copyright © 2000 - 2014 Fluke Biomedical

## Test Record

**TEST PASSED**

### Test Performed

Date: 26-06-2024  
 Record: 801055 2024-06-26.mtr  
 Template: IEC 60601-1 - CLI.mtt  
 Template Version: 1.0.0

### Ansur Components Used

Ansur Version 3.1.3  
 Plug-In: ESA620 Version 1.2.6

## Test Setup

### Selections

Service Events Performed	Standards Performed
	IEC 60601

### Device under test

Serial Number	801055	Type	.
Appliance Code	.	Model	.
Group	.	Location	.
Status	.	Address 1	.
Manufacturer	.	Address 2	..

### MTI Data

Test Instrument	Serial Number	Firmware Version
ESA 620	2863062	v2.12

## Signatures

# Test Result

Test Element	Test Type					Fail
IEC 60601-1 - CL1	Auto Sequence					
<p><b>Procedure:</b></p> <p>(1) Connect the DUT to the ESA620 as indicated in the operators manual.</p> <p>(2) Ensure that DUT power is On.</p> <p>(3) Click module setup and specify the patient leads that are to be tested.</p> <p>(4) Connect patient leads as indicated to the right.</p> <p>(5) Click <b>Start Test</b> to perform the safety test.</p>						
<b>Applied Part setup</b>						
#	Applied Part info		Class	Leads		
1	A.P. Code Serial No. Type	(new)	CF	10		
Mains Voltage	Mains Voltage					
Live to Neutral	Mains Voltage Live to Neutral					
<b>Result:</b> Live to Neutral	<b>Value</b> 121	<b>Unit</b> V	<b>High Limit</b>	<b>Low Limit</b>	<b>Standard</b> IEC 60601	
Neutral to Earth	Mains Voltage Neutral to Earth					
<b>Result:</b> Neutral to Earth	<b>Value</b> 0,2	<b>Unit</b> V	<b>High Limit</b>	<b>Low Limit</b>	<b>Standard</b> IEC 60601	
Live to Earth	Mains Voltage Live to Earth					
<b>Result:</b> Live to Earth	<b>Value</b> 121,1	<b>Unit</b> V	<b>High Limit</b>	<b>Low Limit</b>	<b>Standard</b> IEC 60601	
Protective Earth Resistance	Earth Resistance					
<b>Configuration:</b> Test Current: High						
<b>Result:</b> PE Resistance1	<b>Value</b> 0,137	<b>Unit</b> Ohm	<b>High Limit</b> 0,2	<b>Low Limit</b>	<b>Standard</b> IEC 60601	
Insulation Resistance	Insulation Resistance					
<b>Configuration:</b> Test Voltage: 500V						
Mains to Protective Earth	Insulation Resistance Mains to Protective Earth					
<b>Result:</b> Mains to Protective Earth	<b>Value</b> 99999	<b>Unit</b> MOhm	<b>High Limit</b>	<b>Low Limit</b> 2	<b>Standard</b> IEC 60601	
Applied Parts to Protective Earth	Insulation Resistance Applied Parts to Protective Earth					
<b>Result:</b> Applied Parts to Protective Earth	<b>Value</b> 100,2	<b>Unit</b> MOhm	<b>High Limit</b>	<b>Low Limit</b>	<b>Standard</b> IEC 60601	
Earth Leakage Current	Earth Leakage Current					

Fluke Biomedical Ansur Test Report

Test Element	Test Type	Fail
<i>Configuration:</i> Unused Applied Parts: Floating		
Normal Condition	<i>Earth Leakage Current</i> <i>Normal Condition</i>	
<i>Result:</i> Normal Condition	<b>Value</b> 21,6 <b>Unit</b> uAAC+DC <b>High Limit</b> 5000 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Neutral	<i>Earth Leakage Current</i> <i>Open Neutral</i>	
<i>Result:</i> Open Neutral	<b>Value</b> 41,2 <b>Unit</b> uAAC+DC <b>High Limit</b> 10000 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Normal Condition, Reversed mains	<i>Earth Leakage Current</i> <i>Normal Condition, Reversed mains</i>	
<i>Result:</i> Normal Condition, Reversed mains	<b>Value</b> 20,5 <b>Unit</b> uAAC+DC <b>High Limit</b> 5000 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Neutral, Reversed Mains	<i>Earth Leakage Current</i> <i>Open Neutral, Reversed Mains</i>	
<i>Result:</i> Open Neutral, Reversed Mains	<b>Value</b> 41,4 <b>Unit</b> uAAC+DC <b>High Limit</b> 10000 <b>Low Limit</b>	<b>Standard</b> IEC 60601
<i>Enclosure Leakage Current</i> <i>Configuration:</i> Unused Applied Parts: Floating		
Normal Condition	<i>Enclosure Leakage Current</i> <i>Normal Condition</i>	
<i>Result:</i> Normal Condition	<b>Value</b> 0,9 <b>Unit</b> uAAC+DC <b>High Limit</b> 100 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Neutral	<i>Enclosure Leakage Current</i> <i>Open Neutral</i>	
<i>Result:</i> Open Neutral	<b>Value</b> 0,8 <b>Unit</b> uAAC+DC <b>High Limit</b> 500 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Earth	<i>Enclosure Leakage Current</i> <i>Open Earth</i>	
<i>Result:</i> Open Earth	<b>Value</b> 21,6 <b>Unit</b> uAAC+DC <b>High Limit</b> 500 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Normal Condition, Reversed mains	<i>Enclosure Leakage Current</i> <i>Normal Condition, Reversed mains</i>	
<i>Result:</i> Normal Condition, Reversed mains	<b>Value</b> 0,8 <b>Unit</b> uAAC+DC <b>High Limit</b> 100 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Neutral, Reversed Mains	<i>Enclosure Leakage Current</i> <i>Open Neutral, Reversed Mains</i>	
<i>Result:</i> Open Neutral, Reversed Mains	<b>Value</b> 0,8 <b>Unit</b> uAAC+DC <b>High Limit</b> 500 <b>Low Limit</b>	<b>Standard</b> IEC 60601
Open Earth, Reversed Mains	<i>Enclosure Leakage Current</i> <i>Open Earth, Reversed Mains</i>	
<i>Result:</i> Open Earth, Reversed Mains	<b>Value</b> 20,5 <b>Unit</b> uAAC+DC <b>High Limit</b> 500 <b>Low Limit</b>	<b>Standard</b> IEC 60601

Fluke Biomedical Ansur Test Report

Test Element	Test Type				Fail
Patient Leakage Current Configuration: Total Leakage: No Unused Applied Parts: Floating	Patient Leakage Current				
Normal Condition Result: (new)	Value 0,3	Unit uAAC+DC	High Limit 10	Low Limit	Standard IEC 60601
Open Neutral Result: (new)	Value 0,3	Unit uAAC+DC	High Limit 50	Low Limit	Standard IEC 60601
Open Earth Result: (new)	Value 4,1	Unit uAAC+DC	High Limit 50	Low Limit	Standard IEC 60601
Normal Condition, Reversed mains Result: (new)	Value 0,3	Unit uAAC+DC	High Limit 10	Low Limit	Standard IEC 60601
Open Neutral, Reversed Mains Result: (new)	Value 0,3	Unit uAAC+DC	High Limit 50	Low Limit	Standard IEC 60601
Open Earth, Reversed Mains Result: (new)	Value 3,9	Unit uAAC+DC	High Limit 50	Low Limit	Standard IEC 60601
Mains on Applied Parts Configuration: Total Leakage: No Unused Applied Parts: Floating	Mains on Applied Parts				
Single Fault Condition Result: (new)	Value 12,7	Unit uA	High Limit 50	Low Limit	Standard IEC 60601
Single Fault Condition, Reversed Mains Result: (new)	Value 12,7	Unit uA	High Limit 50	Low Limit	Standard IEC 60601
Patient Auxiliary Current Configuration: Unused Applied Parts: Floating	Patient Auxiliary Current				
Normal Condition Result: (new)	Value 0,4	Unit uAAC+DC	High Limit 10	Low Limit	Standard IEC 60601
Open Neutral	Patient Auxiliary Current Open Neutral				

Fluke Biomedical Ansur Test Report

Test Element	Test Type				Fail
Result:	Value	Unit	High Limit	Low Limit	Standard
(new)	0,4	uAAC+DC	50		IEC 60601
<i>Patient Auxiliary Current</i>					
Open Earth	<i>Open Earth</i>				
Result:	Value	Unit	High Limit	Low Limit	Standard
(new)	2,9	uAAC+DC	50		IEC 60601
<i>Patient Auxiliary Current</i>					
Normal Condition, Reversed mains	<i>Normal Condition, Reversed mains</i>				
Result:	Value	Unit	High Limit	Low Limit	Standard
(new)	0,4	uAAC+DC	10		IEC 60601
<i>Patient Auxiliary Current</i>					
Open Neutral, Reversed Mains	<i>Open Neutral, Reversed Mains</i>				
Result:	Value	Unit	High Limit	Low Limit	Standard
(new)	0,4	uAAC+DC	50		IEC 60601
<i>Patient Auxiliary Current</i>					
Open Earth, Reversed Mains	<i>Open Earth, Reversed Mains</i>				
Result:	Value	Unit	High Limit	Low Limit	Standard
(new)	2,7	uAAC+DC	50		IEC 60601